Detection Module for Chemiluminescence

Jess, Abby, Wes, Peggy Sue and Sally Sue DM-001-DM-006





Let's get started!

These Detection Modules provide all the reagents needed to detect your target of interest with a primary antibody in a Simple Western size-based assay when paired with a Separation Module.

Components

CLAMSHELL BOX

INCLUDES	PART NO
Luminol-S (2 mL)	043-311
Peroxide (2 mL)	043-379
Streptavidin-HRP (132 μL)	042-414
Antibody Diluent 2 or Milk-Free Antibody Diluent (20 mL)	042-203 or 043-524
And one of the following*: Gt Anti-Rabbit Secondary HRP Antibody (2 mL) (min X: n/a)	DM-001
Gt Anti-Mouse Secondary HRP Antibody (2 mL) (min X: n/a)	DM-002
Secondary Streptavidin HRP (2 x 1 mL) ¹ (min X: n/a)	DM-004
Gt Anti-Human IgG Secondary HRP Antibody (2 x 1 mL) ²	DM-005
(min X: HulgA, HulgM, Bov, Hrs, Ms, Rb) Dk Anti-Goat Secondary HRP Antibody (2 x 1 mL) ³ (min X: Chk, GPig, Hrs, Hu, Rb, Rat, Ms, SyrHms)	DM-006

^{*}No secondary antibody is provided with the "No Secondary Detection Module" (PN DM-003).

Storage conditions

Store clamshell at 2-8 °C upon receipt.

Other things you'll need

- Reagents and consumables supplied in the Separation Module.
- An insert with instructions for preparing reagents and pipetting your plate
 is provided with the Separation Module, or it can be downloaded from our
 website.

 See reverse →

¹Use PN 043-459-2 to reorder the Secondary Streptavidin HRP a la carte.

²Use PN 043-491-2 to reorder the Anti-Human IgG Secondary HRP antibody a la carte.

³Use PN 043-522-2 to reorder the Anti-Goat Secondary HRP antibody a la carte.

Detection Module preparation

A ANTIBODIES

- Primary Antibody: Dilute in Antibody Diluent 2. Refer to training material to determine optimal dilution. If you're using the System Control, add 25 µL of the System Control Primary Antibody to 225 µL of your diluted Primary Antibody.
 - If using Goat Primary Antibody, dilute with Milk-Free Antibody Diluent.
- Secondary HRP Conjugate (if applicable): It is provided ready to use in the Detection Module. If you're using your own, consult your FAS or the Antibody Database on the ProteinSimple website.

B SUBSTRATE

• Combine 200 μL Luminol-S and 200 μL Peroxide supplied in the detection module in a microcentrifuge tube.



See the Separation Module product insert for instructions to complete experiment setup.

