

## SIMPLE WESTERN CERTIFIED ANTIBODY DATASHEET

[View  
Antibody Link](#)



Figure1 Simple Western: AIF Antibody [NBP2-27177] - Simple Western lane view shows a specific band for AIF in 0.2 mg/ml of Jurkat and MCF-7 lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system. \*Non-specific interaction with the 230 kDa standard may be seen with this antibody.

PROTEIN TARGET/ANTIBODY	
Protein Target	AIF
Protein Isoform	AIF
Antibody Type	Primary
Host Species/Clonality	Rabbit Polyclonal
ASSAY	
Sample Type	Jurkat and MCF-7 lysate 0.2 mg/ml
Sample Concentration	0.2 mg/ml
Antibody Concentration/Dilution	1:500
Antibody Diluent	
Detection Mode	Chemiluminescence
Separation Type	Size
Matrix	12-230 kDa
Observed kDa	
Publications	None listed

This antibody is validated for Gel-Free, Blot Free, Hands Free Simple Western Systems. To learn about Simple Western Systems, available Simple Western antibodies, or new antibody submissions visit the links below. For additional information, please contact [support@proteinsimple.com](mailto:support@proteinsimple.com).

[Simple Western Systems](#)

[Simple Western Antibody Database](#)

[Simple Western Antibody Submission](#)