


## Simple Western Certified Antibody Datasheet



Simple Western™ analysis of extracts (1 mg/mL) from Human Colon tissue using α-Smooth Muscle Actin (D4K9N) XP® Rabbit mAb #19245. The virtual lane view (left) shows a single target band (as indicated) at 1:10 and 1:50 dilutions of primary antibody. The corresponding electropherogram view (right) plots chemiluminescence by molecular weight along the capillary at 1:10 (blue line) and 1:50 (green line) dilutions of primary antibody. This experiment was performed under reducing conditions on the Jess™ Simple Western instrument from ProteinSimple, a Bio-Techne brand, using the 12–230 kDa Separation Module.

Certified by Cell Signaling Technology

### Antibody

Name	α-Smooth Muscle Actin (D4K9N) XP® Rabbit mAb
Host	Rabbit
Clonality	Monoclonal
Type	Primary
Dilution	1:10–1:50 (recommended by CST)
Sample Concentration	1 mg/mL (recommended by CST)

### Assay

Type	Size
Matrix	12–230 kDa
Observed kDa	48, 49
Sample Type	Human Colon
Publications (2 found)	Nephron. 2022;146(1):110-120. // Int J Mol Sci. 2022 May 9;23(9):5261.

For additional information on this antibody [view antibody link](#).

This antibody is certified for Simple Western™ technology. To learn about Simple Western technology, available antibodies, or to submit new antibodies, visit the links below. For additional information, please contact:

[support@proteinsimple.com](mailto:support@proteinsimple.com)

[Simple Western Systems](#)

[Simple Western Antibody Database](#)

[Submit Antibody Validation Data](#)